

# CAPP® Rondo

## GEL CARD CENTRIFUGE

The CAPP Rondo Gel Card Centrifuge is designed for universal compatibility, making it suitable for use with all standard gel cards. It is an very good choice for hospitals, blood banks and laboratories.

### Features

- This excellent centrifuge is capable of holding up to 12 gel cards simultaneously and is compatible with a wide variety of Gel Cards
- Favourite Mode allows users to set up and store up to 99 user-defined programs
- For enhanced safety, the CAPP Rondo Gel Card Centrifuge features a lid lock mechanism that secures the lid during operation
- Safe operation thanks to automatic imbalance detection
- Microprocessor-controlled centrifuge with variable speed and time setting
- The brushless DC motor enables maintenance-free and long-lasting operation
- Compact design with low loading height for easy operation
- Precise RPM setting from 500 to 2200 RPM in steps of 10 RPM
- Safety lid brake, preventing the lid from opening until the rotor has completely stopped
- Functions - last run memory / RPM/RCF speed setting / short spin



**2 YEAR**  
Warranty

### Specifications

<b>Speed</b>	500-2200 RPM / RCF varies	<b>Noise Level</b>	< 65 dB (A)
<b>Speed Accuracy</b>	± 100 RPM	<b>Time Setting</b>	1 to 999 mins and Infinite mode
<b>RCF</b>	475g		
<b>Motor Type</b>	Brushless DC Motor	<b>Display</b>	Yes
<b>Rotor Capacity</b>	12 Gel cards	<b>Dimension (W x D x H)</b>	381 x 480 x 185 mm
<b>Acceleration Time</b>	35 ± 5 seconds	<b>Weight</b>	8.3kg
<b>Deceleration Time</b>	50 ± 5 seconds	<b>Power Consumption</b>	48 W

## Ordering Information

### Cat. No. Description

CRC-212GC CAPP Rondo Gel Card Centrifuge 2200 RPM

### Scope of Delivery

1 pc. Gel Card Centrifuge

1 pc. Power Adapter

1 pc. T-Allen Key

Reduction Gel Card Plate

Gel Card Holder

1 pc. User Manual

1 pc. Warranty Card