

CAPP Rondo

PLATE TUBE ROTATOR

CAPP Rondo Plate Tube Rotator는 뛰어난 성능을 제공하며 편의성과 신뢰성을 모두 보장합니다. 조용하게 작동하는 브러시리스 DC 모터가 장착되어 있어 유지보수가 필요 없으며 내구성이 매우 뛰어납니다.

Features

- 마이크로프로세서 제어 방식으로, CAPP Rondo Plate Rotator는 가변 속도 설정을 지원합니다.
- 사용자는 최대 9개의 맞춤형 프로그램을 설정하고 저장하여 효율적으로 사용할 수 있습니다.
- 3가지 플레이트 옵션을 지원하며 50mL 튜브 16개, 15mL 튜브 35개, 2mL 튜브 55개를 수용할 수 있어 다양한 실험 용도에 유연하게 대응합니다.
- 자동 전원 복구 기능(Auto power-on)이 있어 정전 후에도 로테이터가 자동으로 다시 작동을 이어서 남은 작업을 완료합니다.



Specifications

Speed	10-80 RPM, in steps of 1 RPM
Speed Accuracy	± 1 RPM
Motor	BLDC motor
Load Capacity	1.5kg
Capacity (supplied with any one chosen disc plate)	16 tubes x 50 mL 35 tubes x 15 mL 55 tubes x 2 mL
Time Setting	1-99 mins and infinite
Pulse Mode	Yes, clockwise and anticlockwise motion with 2 sec stops
Protection class	IP 21
Angle of Disc Plate	Adjustable from 0-90 degrees
Digital Display	Yes
Dimensions (W x D x H)	540 x 274 x 186 mm
Weight	4.5kg
Power Consumption	4.5 W

Scope of Delivery

1 pc. Plate Tube Rotator
1 pc. Power Adaptor
1 pc. User Manual
1 pc. Warranty Card

Ordering Information

Cat. No.	Description
CRR-0800P	CAPPRondo Plate Tube Rotator 80 RPM (plate rotator must be ordered separately)
CRR-0800P-P50	50 mL plate rotator (with clamps)
CRR-0800P-P15	15 mL plate rotator (with clamps)
CRR-0800P-P02	2 mL plate rotator (with clamps)
CRR-0800P-C50	50 mL clamp
CRR-0800P-C15	15 mL clamp
CRR-0800P-C02	2 mL clamp
CRR-0800P-FN	Front nut
CRR-0800P-PA	Power adaptor