

RANDOX
TOXICOLOGY

evidence
MULTISTAT




Rapid Automatic Drug Testing

Evidence MultiSTAT

Using our revolutionary Biochip Array Technology, the Evidence MultiSTAT is an automated analyser that enables on-site simultaneous detection of up to 21 classical, prescription and synthetic drugs from a single sample. Designed to work across a variety of matrices, our patented multi-analyte testing platform provides a complete immunoassay profile within minutes, changing the landscape of drug detection forever.

Benefits



Biochip Array Technology

A patented multi-analyte testing platform provides a complete immunoassay profile.



Multiple matrices

Applicable to a variety of matrices including urine, blood and oral fluid.



Small sample volume

Only a small sample volume required, leaving more for confirmatory testing.



Versatile

Ideal for workplace, custodial and formal laboratory drug testing - a simple interface and pre-filled reagent cartridges allows staff to operate the analyser in any environment.



Automated process

1. Collect and prepare sample
2. Add sample to cartridge
3. Load sample and tip cartridges
4. Press start to get results



Results

Simultaneously screen up to 21 classical, prescription and synthetic drugs, including 4 New Psychoactive Substances, in under 20 minutes from a single sample.

Urine Panel

Cat No EV4193

Assay		
α-PVP (Flakka)	Buprenorphine	Opiates
AB-PINACA	Cannabinoids (THC)	Oxycodone
Amphetamine	Creatinine	Tramadol
Barbiturates	EtG	Tricyclic Antidepressants (TCA)
Benzodiazepines I (Oxazepam)	Fentanyl	URI44 / XLR11
Benzodiazepines II (Lorazepam)	Methadone	JWH-018
Benzoyllecgonine (Cocaine Metabolite)	Methamphetamine	6-MAM

Whole Blood Panel

Cat No EV4195

Assay		
α-PVP (Flakka)	Buprenorphine	Opiates
AB-PINACA	Cannabinoids (THC)	Oxycodone
Amphetamine	Dextromethorphan	Tramadol
Barbiturates	EtG	Tricyclic Antidepressants (TCA)
Benzodiazepines I (Oxazepam)	Fentanyl	URI44 / XLR11
Benzodiazepines II (Lorazepam)	Methadone	JWH-018
Benzoyllecgonine (Cocaine Metabolite)	Methamphetamine	6-MAM

Oral Fluid Panel*

Cat No EV4197

Assay		
α-PVP (Flakka)	Cannabinoids (THC)	Oxycodone
Amphetamine	Fentanyl	Phencyclidine (PCP)
Barbiturates	Ketamine	Tramadol
Benzodiazepines I (Oxazepam)	LSD	URI44 / XLR11
Benzodiazepines II (Lorazepam)	Methadone	JWH-018
Benzoyllecgonine (Cocaine Metabolite)	Methamphetamine	6-MAM
Buprenorphine	Opiates	-

*Optimised for a range of Oral Fluid collection devices

For forensic use only

Drug testing has changed forever